

Most Advanced Access System

- ▶ Advanced seal design
- ▶ Smooth instrument exchange
- ▶ Innovative switch design

Kii[®] Shielded Bladed

Access System



Kii® Shielded Bladed Access System

Shield On



Shield Off



Innovative Switch Design

Provides clear indication of shield position

Funneled Entry

Ensures smooth instrument exchanges with minimal drag



Advanced Seal Design

Offers high flow insufflation and exceptional durability

Facilitates simple seal detachment for rapid desufflation

Optimised Cannula

Accommodates over-sized instruments

Maintains superior cannula retention without barbs

Retention disk allows fixation of the trocar at specific depths



REORDER NO.		DESCRIPTION	QUANTITY
ADVANCED FIXATION	Z-THREAD		
CFB03	CTB03	5mm x 100mm Kii Shielded Bladed Access System	6/BOX
CFB33	CTB33*	11mm x 100mm Kii Shielded Bladed Access System	6/BOX
CFB73	CTB73*	12mm x 100mm Kii Shielded Bladed Access System	6/BOX

REORDER NO.		DESCRIPTION	QUANTITY
NON-THREADED	THREADED		
C0Q60	C0Q61	8mm x 100mm Kii Shielded Bladed Access System	6/BOX

**Please contact your Sales Representative for availability and additional configurations*

Applied Medical
A New Generation Medical Device Company

Rancho Santa Margarita, CA 92688
Tel 949.713.8300
Fax 949.713.8400
www.appliedmedical.com

To place an order, or for more information, call 800.282.2212



© 2011 Applied Medical Resources Corporation. All rights reserved.
Applied Medical, the Applied Medical logo design, ProTerra, the ProTerra logo design and Kii are trademarks of Applied Medical Resources Corporation, registered in one or more of the following countries: Australia, Canada, Japan, the United States and/or the European Community. SC01348J

