

The Last Word in First Entry

- ▶ Establishes pneumoperitoneum under direct visualization
- ▶ Least penetration necessary to initiate insufflation





Kii[®] Fios[®]


First Entry





Kii® Fios® First Entry

Fios First Entry	Insufflation Needle	5mm Access System	12mm Access System
 <ul style="list-style-type: none"> • Least penetration to initiate insufflation • Concurrent direct visualization and rapid insufflation capability 			



Insufflation aperture

Separator Obturator

Enables true non-bladed access with minimal fascial defect



Initial withdrawal



Midway



Complete withdrawal

REORDER NO.		DESCRIPTION	QUANTITY
ADVANCED FIXATION	Z-THREAD		
CFF05	N/A	5 x 75mm Kii Fios First Entry	6/BOX
CFF03	CTF03	5 x 100mm Kii Fios First Entry	6/BOX
CFF33	CTF33	11 x 100mm Kii Fios First Entry	6/BOX
CFF73	CTF73	12 x 100mm Kii Fios First Entry	6/BOX

Available in a full range of sizes and configurations.



Rancho Santa Margarita, CA 92688
 Tel 949.713.8300
 Fax 949.713.8400
www.appliedmedical.com

To place an order, or for more information, call 800.282.2212



© 2011 Applied Medical Resources Corporation. All rights reserved.
 Applied Medical, the Applied Medical logo design, ProTerra, the ProTerra logo design and marks designated with a ® are trademarks of Applied Medical Resources Corporation, registered in one or more of the following countries: Australia, Canada, Japan, the United States and/or the European Union. SC01383H

