### One Hand. One Step.

- ► Easy duct access
- ▶ Self-anchoring
- ► One-handed activation

# <u>AeroStat®</u>

Cholangiography Catheter







## AeroStat® Cholangiography Catheter

#### **Easy Duct Access**

Retention cone eliminates problems traditionally associated with laparoscopic cholangiography, such as clipping, balloon purging and potential breakage.

#### Self-anchoring

Soft 6F-20F expandable cone and 3F tapered rigid tip atraumatically engages and occludes the cystic duct, eliminating the need for clips or graspers.

#### **One-handed Activation**

Ratcheting handle, with one-handed activation, makes placement fast and easy.

REORDER NO.	DESCRIPTION	SHAFT DIAMETER	LENGTH
C1002	Rigid Tip with Percutaneous Introducer	6F	33cm



Rancho Santa Margarita, CA 92688 Tel 949.713.8300 Fax 949.713.8400

To place an order, or for more information, call 800.282.2212

www.appliedmedical.com



