


# Omni-Flex™

General Surgery Retractor System



Your Patient.  
Your Approach.

## Flexibility

### Interchangeable, Removable Arms

Stable retraction that lets you perform procedures the way you prefer. Strong enough for your largest patients. Nimble enough for small-site work. Always easy to use.

Use one arm. Use two arms. Take one arm away or use one of each. Choose the frame arm—before or during surgery—that helps maximize your surgical field.

## Ease-of-Use

### Maximize Exposure—Minimize Obstructions

Exchange an arm. Remove an arm. Open the field. The Omni-Flex System easily exposes difficult anatomy no matter what your procedure demands.

## Strength

### Strength for Stability with Every Patient

Any patient size. Any incision size. Always stable. The Omni-Flex System is designed to remain stable, so you always get the exposure you need.





Retract from every direction by simply adding our hinged extension arm.



*flexibility*

*ease-of-use*



**Flexibility for**

Any approach. Any plane. Any position. Always low-profile. The O



## SideSnap™ Clamps for Faster Setup & Mid-Surgery Changes

Snap the clamp on the frame. Snap the blade in place. Our SideSnap™ Clamps allow for quicker, simpler setup. With added strength, they enhance the flexibility of the Omni-Flex System.



## Quick Blade Attachment

Attach in one sweeping motion with a simple snap from the side.



## Stack Blades to Maximize Incision

Stack a longer blade over a shorter blade for increased tissue holding power and increased visibility.




of-use  
strength



## Any Approach

Omni-Flex Retractor System gives you complete surgical freedom.

# Product Specifications



## Omni-Flex™ General Retractor Systems

### G0550 Expanded General Retractor System

### G0500 General Retractor System

G0550	G0500	Description	Part Number	
A	1	1	Sterile Field Post	10244
B	1	1	Support Arm	11050
C	2	2	Standard WISHBONE® Frame Arm	10269
D	2	2	Small WISHBONE® Frame Arm	10270
E	8	6	SideSnap™ Clamp	11030
F	1	--	Hinged Extension Arm, 19" (48.3cm)	11121
G	1	1	Deaver Swivel Retractor, 2 3/4" x 4 1/2" (7cm x 11.4cm)	3784
H	1	1	Harrington Sweetheart Swivel Retractor, 2" x 5" (5cm x 12.7cm)	3794
I	1	1	Fence Swivel Retractor, 4" x 6" (10.2cm x 15.2cm)	3772
J	1	1	Kelly Modified Swivel Retractor, 3" x 2 1/2" (7.6cm x 6.4cm)	3862
K	1	1	Kelly Modified Swivel Retractor, 3" x 4" (7.6cm x 10.2cm)	3864
L	1	1	Malleable Slotted Swivel Retractor, 2" x 6 1/2" (5cm x 16.5cm)	3788
M	1	1	Malleable Slotted Swivel Retractor, 3" x 6 1/2" (7.6cm x 16.5cm)	3790
N	2	2	Mayo Swivel Retractor, 2 3/4" x 2" (7cm x 5cm)	3335
O	1	1	Splanchnic Swivel Retractor, 2" x 6" (5cm x 15.2cm)	3774



4849 White Bear Parkway, St. Paul, MN 55110 USA  
 800-367-8657 • 651-287-4300 • fax: 651-287-4400  
[www.omni-tract.com](http://www.omni-tract.com)

© 2006 Omni-Tract Surgical  
 All Rights Reserved.



3-06/40-6154/1000