



CONFIDENCE. SECURITY. HEM-O-LOK.®

WECK® Hem-o-lok polymer locking ligation system.

AT A GLANCE

- Provides secure ligation, clip locks on to vessels
- Tactile feedback
- Permanent, non-absorbable, inert, non-toxic, radiolucent polymer
- Cool ligation system
- Open or laparoscopic appliers

CONFIDENCE YOU CAN FEEL

Confidence. That's what the WECK Hem-o-lok polymer locking ligation system provides. Confidence created by a ligation system designed for clip security—secure in the applier, locked on to the patient vessel.

SECURITY BY DESIGN

Hem-o-lok's unique design and flexible hinge keep the clip firmly seated in applier jaws. The clip's distinct tactile feedback confirms jaw seating and secure vessel placement.

COOL LIGATION


Since Hem-o-lok is a cool ligation system, there's no chance of thermal spread to vital structures, a common concern of many surgeons when they use energy-based ligation solutions.

SIZE RANGES TO MEET YOUR NEEDS

Hem-o-lok is available in a range of clip sizes (Medium, Medium/Large, Large, and XL) to ligate up to 10 mm vessels through a 5 mm port or up to 16 mm vessels through a 10 mm port. Our broad range of applier lengths (including 45 cm endo appliers for extended reach) allows surgeons to deliver clips with control and confidence.



VESSEL/TISSUE BUNDLE SIZE RANGE




HEM-O-LOK POLYMER LIGATING CLIPS

- 544220 Medium
- 544230 Medium-Large
- 544240 Large
- 544250 Extra Large

CLIPS	OPEN APPLIERS/REMOVERS	ENDOSCOPIC APPLIERS/ REMOVERS	AUTOMATIC APPLIERS
 <p>M</p> <p>544220 Hem-o-lok® Medium Polymer Clips 6 clips/cartridge 14 cartridges per box, sterile</p>	<p>544113 8" Curved 544114 11" Right-Angle (70°) 544115 11" Curved 544123 11" Remover (for M, ML, L)</p>		
 <p>ML</p> <p>544230 Hem-o-lok® Medium-Large Polymer Clips 6 clips/cartridge 14 cartridges per box, sterile</p>	<p>544170 8" Curved 544171 11" Curved 544172 11" Right-Angle (70°) 544123 11" Remover (for M, ML, L)</p>	<p>Endo5® ML, 5mm diameter 544965 32.5cm applier 544965L 45cm applier 544965T 34cm applier, take-apart 544121 32cm remover (for ML and L) 544121L 45cm remover (for ML and L) 544121T 32cm remover, take-apart (for ML and L)</p>	<p>543965 Auto Endo5® 15 clips per applier, 3 appliers per box, sterile, disposable for endo surgery</p>
 <p>L</p> <p>544240 Hem-o-lok® Large Polymer Clips 6 clips/cartridge 14 cartridges per box, sterile</p>	<p>544179 11" Right-Angle (70°) 544180 8" Curved 544181 11" Curved 544123 11" Remover (for M, ML, L) 544124 11" Remover (for L, XL)</p>	<p>Endo10™ML 544995 32cm applier, 10mm 544995L 45cm applier, 10mm 544995T 34cm applier, 10mm, take-apart 544121 32cm remover, 5mm (for ML and L) 544121L 45cm remover, 5mm (for ML and L) 544121T 32cm remover, 5mm, take-apart (for ML and L) 544130 32cm remover, 10mm (for L and XL) 544130L 45cm remover, 10mm (for L and XL) 544130T 33cm remover, 10mm take-apart (for L and XL)</p>	<p>543181 Auto Open 15 clips per applier, 3 appliers per box, sterile, disposable for open surgery</p>
 <p>XL</p> <p>544250 Hem-o-lok® Extra-Large Polymer Clips 6 clips/cartridge 14 cartridges per box, sterile</p>	<p>544191 11" Curved 544192 11" Right-Angle (70°) 544124 11" Remover (for L, XL)</p>	<p>Endo10™ XL, 10mm diameter 544990 32cm applier 544990L 45cm applier 544990T 34cm applier, take-apart 544130 32cm remover (for L and XL) 544130L 45cm remover (for L and XL) 544130T 33cm remover, take-apart (for L and XL)</p>	
STERILITY STORAGE TRAY		544500 capacity: 3 endoscopic 32cm appliers	P37319 capacity: 2-3 endoscopic 45cm appliers

HEM-O-LOK BEHIND THE CLIP



- Bosses** are designed to retain clip in applier jaws
- Integrated teeth** interface with the vessel and are designed to prevent slippage
- Locking mechanism** provides tactile feedback and secure closure
- Bow-shape** allows removal with appropriate instrument
- Hinge** allows flexibility in clip placement prior to clip locking
- Non-Absorbable polymer** is inert, non-conductive, radiolucent, and does not interfere with CT, MRI, or X-ray diagnostics

TO PLACE AN ORDER CONTACT CUSTOMER SERVICE AT 1.866.246.6990

TELEFLEX MEDICAL Research Triangle Park, NC 27709
Phone: 866.246.6990 Fax: 800.399.1028 Intl: 919.361.3961
WWW.TELEFLEXMEDICAL.COM

Hem-o-lok, Weck, Endo5 and Teleflex are registered trademarks of Teleflex Incorporated or its affiliates.
©2010 Teleflex Incorporated. All rights reserved. 2010-0050 v1

