

ORDERING INFORMATION

□ Micro Small ■ Small/Wide ■ Medium ■ Medium/Large ■ Large

TITANIUM LIGATING CLIPS

Product Number	Clip Size	Clips/ Cartridge	Cartridges/ Box	Clips/ Box
005200	Microclip	6	30	180
001200	Small	6	30	180
001204	Small	24	25	600
001201	Small-Wide	6	30	180
001205	Small-Wide	24	25	600
002200	Medium	6	30	180
002204	Medium	24	25	600
003200	Medium-Large	6	20	120
003204	Medium-Large	24	15	360
004200	Large	6	20	120
004204	Large	24	15	360

ENDOSCOPIC LIGATING APPLIERS

Product Number	Applier Length	Diameter	Clip Size	Jaw Type
238110	11" (28cm)	10mm	Medium	Endoscopic
238170	18" (45cm)	10mm	Medium	Endoscopic
238200	13" (32cm)	5mm Shaft, 7mm Jaw	Medium	Endosaphenous
338110	11" (28cm)	10mm	Medium-Large	Endoscopic
338170	18" (45cm)	10mm	Medium-Large	Endoscopic
438110	11" (28cm)	12mm	Large	Endoscopic
438170	18" (45cm)	12mm	Large	Endoscopic

SPECIALTY OPEN SURGICAL LIGATING APPLIERS

137117	11" (28cm)	Small-Wide	Multi-Angled	45° Angled
137118	11" (28cm)	Small-Wide	Multi-Angled	Curved

METAL CLIP REMOVERS (REMOVES ALL WECK METAL CLIPS)

Product Number	Applier Length	Style
523120	6.5" (16.5cm)	Straight
523121	10.5" (26.7cm)	Bent Shaft, 30° Angle

OPEN SURGICAL LIGATING APPLIERS

Product Number	Applier Length	Clip Size	Jaw Type*
537061	6" (15cm)	Microclip	Curved
537081	8" (20cm)	Microclip	Curved
137061	6" (15cm)	Small	Curved
137081	8" (20cm)	Small	Curved
137085	8" (20cm)	Small	Angled
137111	11" (28cm)	Small	Curved
137062	6" (15cm)	Small-Wide	Curved
137082	8" (20cm)	Small-Wide	Curved
137086	8" (20cm)	Small-Wide	Angled
137112	11" (28cm)	Small-Wide	Curved
237061	6" (15cm)	Medium	Curved
237081	8" (20cm)	Medium	Curved
237085	8" (20cm)	Medium	Angled
237111	11" (28cm)	Medium	Curved
237115	11" (28cm)	Medium	Angled
237117	11" (28cm)	Medium	Right-Angled
237141	14" (35.5cm)	Medium	Curved
237145	14" (35.5cm)	Medium	Angled
337081	8" (20cm)	Medium-Large	Curved
337085	8" (20cm)	Medium-Large	Angled
337111	11" (28cm)	Medium-Large	Curved
337115	11" (28cm)	Medium-Large	Angled
337117	11" (28cm)	Medium-Large	Right-Angled
337141	14" (35.5cm)	Medium-Large	Curved
337145	14" (35.5cm)	Medium-Large	Angled
437081	8" (20cm)	Large	Curved
437085	8" (20cm)	Large	Angled
437111	11" (28cm)	Large	Curved
437115	11" (28cm)	Large	Angled
437117	11" (28cm)	Large	Right-Angled
437141	14" (35.5cm)	Large	Curved
437145	14" (35.5cm)	Large	Angled

*Curved 20° - Angled 55° - Right-Angle 70°

WECK® HORIZON™ LIGATION SYSTEM

The Name You Know
The Product You Trust



TO PLACE AN ORDER CONTACT CUSTOMER SERVICE AT 1.866.246.6990

TELEFLEX MEDICAL Research Triangle Park, NC 27709
Phone: 866.246.6990 Fax: 800.399.1028 Intl: 919.361.3961
WWW.TELEFLEXMEDICAL.COM

All Horizon™ ligating clips are packaged sterile. Ligating clip systems differ in their clip-closure characteristics. It is the responsibility of the user to confirm security of clips on vessels or other tissue. Not intended for contraceptive tubal occlusion.

Weck® is a registered trademark and Horizon™ and Microclip™ are trademarks of Teleflex Incorporated or its affiliates. ©2010 Teleflex Incorporated. All rights reserved. 2008-0085 Rev 01

Microclip™
Horizon's newest
family member.



For safe ligation of small
vessels or structures.



Cross-sectional view of
clip with applicator jaws
properly aligned (above)
and misaligned (below).



THE HORIZON SYSTEM: DESIGNED TO MEET YOUR NEEDS



Horizon is designed to meet the needs of surgeons, nurses and other operating room personnel.

THE CLIP: FUNCTIONAL EXCELLENCE

- The chevron shape of the Horizon clip engulfs tissue with precise tip-to-tip closure.
- Heart-shaped wire is designed to give each clip a firm grip on vessels.
- Transverse grooves make the clip resistant to slipping across the vessel.
- The triangulated cross-section of the Horizon clip leg maximizes surface-to-surface contact of clip and applicator jaw, practically eliminating clip fallout.
- The internal groove of the Horizon clip is sized to reduce clip scissoring.



THE CARTRIDGE: EASY LOADS WITH ADHESIVE CONVENIENCE

- The Horizon clip cartridge can provide incredibly easy and precise loading that can enhance surgical procedure efficiency.
- Lateral clip-restraining springs in the Horizon cartridge are displaced by applicator jaws during clip loading.
- Spring displacement frees the clip for a nearly frictionless one-hand load from the cartridge.
- Clips are centered in the cartridge slots for precise loads.



- Clip cartridges are color-coded to match applicators.
- Both 6- and 24-clip cartridges are available.
- Special adhesive backing on the Horizon cartridge allows the user to conveniently relocate the site of clip loading — cartridges affix to any sterile area.



THE APPLIERS: PERFORMANCE AND CLEANING ACCESSIBILITY

- The Horizon endoscopic applicator provides clip closure designed to give surgeons solid tactile feedback.
- The large, lightweight rotation knob allows for one-finger, 360° rotation of the applicator shaft.
- Horizon endoscopic applicators are available for Medium, Medium-Large and Large Horizon clips and feature flush ports for thorough cleaning.
- The reusable applicators are designed specifically for use with Horizon ligating clips, creating a dramatic cost savings when compared to disposable automatic-feed clip applicators.
- All Horizon open surgical ring-handled applicators open at the box lock for easy, thorough cleaning.
- Open surgical applicators are available with angled, curved or right-angled jaws.



Specialty Multi-Angle Applicators

The Horizon 11" Multi-Angle Access Applicator for small (wide slot) clip cartridges incorporates 45° angled jaw, or a curved jaw, for an optimal angle of access to the surgical target.



Specialty Endoscopic Applicators

The Horizon ESV Applicator incorporates a 5mm shaft diameter (7mm jaw) with an angled jaw for ligation of vein branches.

